


supplement

**EMI**

# photomultiplier tubes





# Photomultiplier Tubes

## Catalogue Supplement

### INDEX

#### List of tube types (including variants)

Type	Page	Type	Page
9558F	15	9813B	25
9633B	11	9813KB	25
9633KB	11	9813QB	25
9633QKB	11	9813QKB	25
9633QKB	11	9814B	25
9661B	5	9814KB	25
9664B	5	9814QB	25
9665B	5	9814QKB	25
9670B	5	9815B	25
9683B	11	9815KB	25
9683KB	11	9815QB	25
9683QB	11	9815QKB	25
9683QKB	11	9816B	25
9757B	17	9816KB	25
9769B	31	9816QB	25
9778B	21	9816QKB	25
9778KB	21	9817B	25
9781B	5	9817KB	25
9781R	5	9817QB	25
9783B	5	9817QKB	25
9783R	5	9818B	25
9785B	5	9818KB	25
9785QB	5	9818QB	25
9789B	19	9818QKB	25
9789QB	19	9820B	23
9804B	17	9820QB	23
9810B	25	9824B	9
9810KB	25	9824QB	9
9810QB	25	9836B	31
9810QKB	25	9840B	31
9811B	25	D125B	9
9811KB	25	D125QB	9
9811QB	25	D214B	13
9811QKB	25	D903B	3
9812B	25		
9812KB	25		
9812QB	25		
9812QKB	25		

#### Sections

	Page
Introductory notes	1
Tubes, diameter less than 2 in	3-12
Side window tubes	5-8
Tubes, diameter 2 in and above	13-29
Fast linear focused tubes	25-29
Magnet assemblies	30
Information on other new tubes and other products	31
Photomultiplier tube power supply P.M.25A	31

#### NOTE

The manufacturer reserves the right to modify designs and specifications without prior notice. Every effort is made to ensure accuracy but no responsibility can be accepted for errors or consequences arising therefrom.

# Introductory Notes

This publication is a supplement to the EMI photomultiplier catalogue reference P001/fP70. It gives data and specifications on the new tube types which have been introduced since the printing of the main catalogue. Some of these tubes are still developmental and their specifications therefore are headed "Provisional Data". Equipment manufacturers are advised to consult EMI before designing such types into production instruments.

The new tubes incorporate many new design features which provide an even wider field of application for the EMI photomultiplier tube range. This is something which should interest all photomultiplier tube users irrespective of their own particular requirements.

## Other EMI publications

The main photomultiplier catalogue P001/fP70 contains a section describing the fundamentals of photomultiplier tube operation and design parameters. Copies are available free on request.

Also available at no charge is a range of technical papers covering a wide variety of specialised subjects. A list is given on the next page.

## Spectral Response Calibration Service

For several years a monochromator spectral calibration service has been available for photomultipliers, denoted by the suffix 'M' following the type number.

New equipment has now been developed which automatically measures the quantum efficiency in per cent and the radiant sensitivity in mA/W at wavelength intervals of 20nm so that the spectral response curve for a particular tube can be plotted in either unit.

The spectral measurements are made by comparing the response of the tube on test with that of a previously calibrated reference photodetector. The reference is calibrated by the same comparison method using a tube whose response has been measured by the National Physical Laboratory.

By integration of the spectral response curve, figures are also derived for the photocathode sensitivity in  $\mu\text{A}/\text{lm}$ , and the Corning Red, Corning Blue and Wratten 87 filter readings where appropriate. These are used by the test engineer to check the standard sensitivity results given on the tube test ticket.

For further information about this service, please contact EMI Electron Tube Division.

## General notes on photomultiplier tube operating conditions

- a) Take great care when clamping tubes, particularly those with quartz windows, as excess pressure may cause fracture and render the warranty void.
- b) Tubes with plain glass bases should be used whenever minimum stray capacitance and electrical leakage are required. Matching Teflon sockets are supplied for this purpose.
- c) Teflon sockets should not be used when the temperature is less than  $-40^\circ\text{C}$ . Individual connectors (Hyperboloidal contacts) are available for this purpose.
- d) Any material in contact with the glass envelope must be maintained at cathode potential. Failure to do so may result in erratic operation and high dark current.
- e) Magnetic shields should be used wherever possible (see page 64 of EMI Photomultiplier Catalogue).
- f) Overall sensitivity is greatly affected by variations in overall voltage, particularly in the case of tubes with a large number of dynodes.
- g) Tubes in the range 9810 to 9820 inclusive incorporate envelope screens which are connected to the photocathode terminal. In many applications the photocathode will be operated at high potential; therefore extreme care should be taken to prevent the possibility of endangering personnel in the vicinity.

For further information the address to contact is EMI Electronics Limited, Electron Tube & Microelectronics Division, 243 Blyth Road, Hayes, Middlesex, England.

## Reprints of Published Technical Papers\*


All reprints are supplied free of charge.

Title	Ref. No.
<b>Neutron Detectors.</b> By J. Sharpe, B.Sc., C.Eng., F.I.E.E.	R/P001
<b>Photomultiplier Tubes and Scintillation Counters.</b> By J. Sharpe, B.Sc., C.Eng., F.I.E.E., and E. E. Thomson, Ph.D., M.Sc., D.I.C.	R/P002
<b>Photoelectric Cells and Photomultipliers.</b> By J. Sharpe, B.Sc., C.Eng., F.I.E.E.	R/P006
<b>Photomultipliers for Tritium Counting. Reprinted from "Tritium in the Physical and Biological Sciences".</b> By J. Sharpe, B.Sc., C.Eng., F.I.E.E., and V. A. Stanley, M.Sc.	R/P007
<b>Vacuum and Gas Photocells and Semiconductor Photocells.</b> By J. Sharpe, B.Sc., C.Eng., F.I.E.E.	R/P009
<b>Fatigue and Saturation in Photomultipliers.</b> By J. P. Keene, B.Sc.	R/P015
<b>A Photomultiplier Circuit for Precision Spectrophotometry.</b> By W. Budde	R/P019
<b>Rapid Ozone Determination near an Accelerator.</b> By G. W. Nederbragt, A. Van Der Horst, and J. Van Duijn	R/P020
<b>Dark Current in Photomultiplier Tubes.</b> By J. Sharpe, B.Sc., C.Eng., F.I.E.E.	R/P021
<b>Photomultiplier Applications to Nucleonic Measurement.</b> By R. B. Owen, B.Sc., M.I.E.E.	R/P024
<b>An Application of Photomultiplier Tubes in Temperature Measurement.</b> By J. Sharpe, B.Sc., C.Eng., F.I.E.E.	R/P025
<b>Method for the automatic measurement of periodically changing spectral energy distributions of extremely weak light intensities.</b> By J. Schanda, Hungary	R/P026
<b>Archaeology—Thermoluminescent Dating of Ancient Ceramics.</b> By M. J. Aitken, M. S. Tite and J. Reid	R/P027
<b>Tests of Photomultipliers for Astronomical Pulse-Counting Applications.</b> By J. P. Rodman and H. J. Smith	R/P028
<b>Measurements of very low Spectral Intensities.</b> By H. A. W. Tothill, Ph.D., A.C.G.I.	R/P029
<b>Photocathodes used in Photomultiplier Tubes for Scintillation Counters—Characteristics and effect on performance.</b> By V. A. Stanley, M.Sc.	R/P031
<b>Some Statistical Properties of Pulses from Photomultipliers.</b> By M. Gadsden	R/P032
<b>EMI Scintillation Counters.</b> By J. Sharpe, B.Sc., C.Eng., F.I.E.E.	R/P033
<b>Notes on the Performance and Application of EMI windowless Particle Multipliers.</b> By A. B. Smith, B.A., A.M.I.E.E.	R/P034
<b>Electron Capture Decay of <sup>40</sup>K.</b> By M. F. McCann, G. M. Lewis and K. M. Smith	R/P037
<b>A synchronized direct reading spectrometer-spark system for time resolved spectra. The study of time-resolved spectra of copper spectra analysis lines used for the analysis of copper in aluminium alloys.</b> By W. W. Schroeder and A. Strasheim	R/P039
<b>Some Characteristics of a Venetian Blind Particle Multiplier used as a Detector in the Study of Ion-Molecule Reactions.</b> By K. R. Ryan	R/P048
<b>Homogeneous Chemiluminescent Measurement of Nitric Oxide with Ozone. Implications for continuous selective monitoring of gaseous air pollutants.</b> By Arthur Fontijn, Alberto J. Sabadell, and Richard J. Ronco, AeroChem Research Laboratories, Inc., P.O. Box 12, Princeton, N.J. 08540—Reprinted from "Analytical Chemistry", Vol. 42, No. 6, May 1970	R/P050
<b>On the use of a Photomultiplier as a Photon Counter.</b> By J. F. James	R/P054
<b>Photomultiplier tubes and their application.</b> By G. C. Young	R/P055



\*Please quote titles and reference numbers when ordering.



OVERALL SENSITIVITY vs OVERALL VOLTAGE



SPECTRAL RESPONSE




All dimensions in millimetres (inches in parenthesis).

**PIN CONNECTIONS** (viewed from below starting left of short pin or key)

Type	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	socket
D903B	—	d1	d3	d5	d7	d9	d11	A	d12	d10	d8	d6	d4	d2	K	flying leads

# Side Window Tubes

- Nominal tube diameter: 30 mm (1 $\frac{1}{8}$  in)
- 9 compact focused dynodes with CsSb secondary emitting surfaces
- Tubes overcapped with base type B11A (B11-88)  
Uncapped versions are available as follows:
  - 9661B available uncapped as 9660B
  - 9665B " " " 9662B
  - 9664B " " " 9663B
  - 9670B " " " 9669B



## PRELIMINARY SELECTION GUIDE

Type	Spectral Response	Wavelength Range (nm)	Window Material	Remarks
9781B	Modified S-5	185-650	UV glass	High performance parent type
9783B	Modified S-5	165-650	Quartz	Spectrosil version of 9781B
9781R	Extended S-5	185-750	UV glass	9781B variant with guaranteed red response
9783R	Extended S-5	165-750	Quartz	Spectrosil version of 9781R
9785B	Trialkali	185-830	UV glass	High performance, wide spectral coverage
9785QB	Trialkali	165-830	Quartz	Spectrosil version of 9785B
9661B	S-5	185-650	UV glass	General purpose, UV detection
9665B	S-5	165-650	Quartz	Quartz version of 9661B
9664B	S-10	185-800	UV glass	I.R. response (see also 9785B)
9670B	S-10	165-800	Quartz	Spectrosil version of 9664B (see also 9785QB)

## RATINGS (ALL TYPES)

### Overall sensitivity:

9781B 9783B 9781R 9783R	Rated at 50 & 200 A/lm
9785B 9785QB	Rated at 1000V
9661B 9665B 9664B 9670B	Rated at 20 & 50 A/lm
Voltage, cathode to d1: maximum	150V
Voltage, anode to cathode: maximum	1250V
Anode current (mean): maximum	0.5mA
Anode dissipation: maximum	0.5W
Anode pulse rise time: typical	2 n. sec.
Anode pulse f.w.h.m.: typical	4 n. sec.
Transit time: typical	20 n. sec.
Capacitance, anode to all dynodes: typical	6pF
Operating temperature: maximum	60°C
minimum	-180°C

## TUBE REPLACEMENT DATA

The following EMI types are direct plug-in replacements, although it may be necessary to adjust the overall voltage to obtain the same operating conditions. If this is the case, it is normally a reduction in overall voltage that is required.

Tube to be replaced		EMI replacement type
RCA	HTV	
1P21	1P21	9781A*
1P28	1P28	9781B, 9661B
1P22		9661B
1P28A		9781R
	R212	9781B
	R106	9783B
	R136	9785QB
	R213	9785B
	R446	9785B
	R456	9785QB

\* 'A' denotes selected tube. When ordering please specify "equivalent to 1P21".

**CHARACTERISTICS**

	Cathode Sensitivity $\mu\text{A}/\text{lm}$		Overall Sensitivity 50A/lm				Overall Sensitivity 200A/lm			
			V		Dark Current nA		V		Dark Current nA	
	Min.	Typ.	Typ.	Max.	Typ.	Max.	Typ.	Max.	Typ.	Max.
9781B 9783B	40	60	670	800	1	10	800	—	4	—
9781R* 9783R*	40	60	670	800	1	10	800	—	4	—
			20A/lm				50A/lm			
9661B 9665B	20	40	750	1000	2	20	850	—	5	—
9664B* 9670B*	20	35	750	1000	2	20	850	—	5	—

\* For red sensitivity specification see adjacent table.

	Cathode Sensitivity $\mu\text{A}/\text{lm}$		Overall Red to Blue Sensitivity Ratio†	Overall Voltage 1000V			
				Overall Sens. A/lm		Dark Current nA	
	Min.	Typ.	Min.	Min.	Typ.	Typ.	Max.
9785B 9785QB	20	50	10%	40	250	5	50

† See adjacent table for method of measurement.

**RED SENSITIVITY**

**9785B, 9785QB**


The red/blue sensitivity ratio is determined by operating the tube at 50A/lm and adjusting the white light intensity to give an anode current of 1mA. A Kodak Wratten 87 filter, which passes all radiation of wavelength longer than 750 nm, is placed between the photocathode and the light source and the resultant anode current measured. The ratio of this reading to a similar reading taken with a Corning blue filter CS-5-58 (polished to half stock thickness) is calculated.

**9781R, 9783R, 9664B, 9670B**




In order to specify the red sensitivity of these types, the tube is operated at 50A/lm and the light intensity adjusted to give an anode current of 1mA. A Kodak Wratten 88A filter, which passes all radiation of wavelength longer than 730nm, is placed between the photocathode and the light source and the resultant anode current in  $\mu\text{A}$  taken as an indication of the red response of the tube.

	Red Sensitivity with Wratten 88A filter	
	Min.	Typ.
9781B	—	0.1
9781R	0.5	2.0
9783R	—	—
9664B	2.0	5.0
9670B	—	—

**SPECTRAL RESPONSE**



OVERALL SENSITIVITY vs OVERALL VOLTAGE



## Notes

### 1 Test Ticket Information

- Each tube is individually calibrated and the test ticket provided with the tube specifies the cathode sensitivity in  $\mu\text{A}/\text{lm}$ , the overall voltages for 50  $\text{A}/\text{lm}$  and 200  $\text{A}/\text{lm}$  (9781, 9783) or 20  $\text{A}/\text{lm}$  and 50  $\text{A}/\text{lm}$  (9661, 9665, 9664, 9670) and the dark currents at these overall sensitivities. A Wratten 88A filter reading is also given for types 9781R, 9783R, 9664B and 9670B (see page 6). For the 9785, the cathode sensitivity is specified in  $\mu\text{A}/\text{lm}$ . The overall sensitivity in  $\text{A}/\text{lm}$  and dark current are given at an overall voltage of 1000V. A red/blue sensitivity ratio is also quoted (see page 6).
- Test data is obtained with cathode—d1 voltage held at 100V and a linear dynode chain. The light source is a tungsten filament lamp operated at a colour temperature of 2857°K. The overall sensitivity is given for illumination of the effective cathode area, which is less than the geometrical area. For the cathode sensitivity measurement, 100V are applied between cathode and d1 with all other electrodes connected to d1.

### 2 Selection

It is possible to select tubes for one or more parameters to meet individual customer's requirements. In this case, the suffix 'A' is used in place of the 'B' or in addition to the 'R' e.g. 9785A, 9781RA.

### 3 Stability

- The magnitude of the anode current determines the operating stability of these tubes and for highest stability in d.c. conditions, the mean anode current should not exceed  $1\mu\text{A}$ . However, mean anode currents up to  $10\mu\text{A}$  will not seriously affect the stability of operation.
- For all tubes, the overall multiplication is extremely sensitive to variations in the applied voltage and the output of the power supply must therefore be very stable and have minimal ripple.
- Any material in contact with the glass envelope must be held at cathode potential. Failure to do so may result in erratic operation and high dark current. (Great care should be taken in clamping tubes, particularly those with quartz windows as excess pressure may fracture the glass and invalidate the guarantee).

### 4 Temperature Considerations


The dark current of all photomultiplier tubes is a very critical function of temperature and operation of tubes above normal ambient temperature ( $22^\circ\text{C}$ ) will result in a higher dark current than that specified on the test ticket. Tubes may be selected for very low dark current but for very stringent applications the use of a cooled photomultiplier housing is recommended.

### 5 Dynode Chain Design

A typical dynode chain is shown below. The dynode resistance value should be sufficient to prevent the dynode potential varying with signal current and a dynode chain current of 10 times the mean anode current is generally suitable. The dynode resistor value, R, is normally chosen to be between  $10\text{K}\Omega$  and  $100\text{K}\Omega$ . The use of a zener diode between k—d1 is not critical and a resistor may be used.

For pulsed applications, when the ratio of peak anode current to mean anode current is high, the last dynode should be decoupled to ground with an  $0.1\mu\text{F}$  capacitor and the next two dynodes decoupled with  $0.01\mu\text{F}$  capacitors.

Further details of dynode chain design and power supply requirements are given in the EMI Photomultiplier Catalogue, reference P001/fP70, available on request.




### 6 Base Configurations


All types are supplied overcapped with a small-shell sub-magnet 11-pin base type B11A (B11-88). A matching PTFE B11A socket is available from EMI and should be installed so that the keyway in the socket faces the incident radiation.

The uncapped tubes, types 9660B, 9662B, 9663B and 9669B are supplied complete with matching PTFE socket type B14B.

#### Tube Types 9781, 9783, 9661, 9665, 9664, 9670, 9785




#### Tube Types 9660, 9662, 9663, 9669,



All dimensions are in millimetres

#### Pin Connections

(Viewed from below)



# New 30mm(1½in) diameter tubes with bialkali or S20 photocathodes and 11 dynodes

## Provisional Data

- Wide spectral coverage —210 to 850 nm (170–850 nm quartz version)
- High gain—typically  $10^6$  at 1000V
- Physically interchangeable with S10 tubes EMI 9592, 9529 and S11, EMI 9524, 9526 and 9601
- Small size—30 mm diameter envelope with 23 mm photocathode

The 9824 (previously D110) is intended to supersede the 9524S in applications requiring a tube of small physical size having a combination of low background and good "blue" response. Such applications include low energy nuclear radiation monitoring, geological studies and low light level photometry.

The D125 is a developmental tube primarily designed as an alternative to existing S10 types where higher red sensitivity is required without circuit changes. Applications include wide range spectrophotometry, photometry, laser detection, colorimetry, film scanning and automated process control.

### Notes

- 1 Each tube is individually calibrated and supplied with a test ticket giving the cathode sensitivity in  $\mu\text{A}/\text{lm}$ ; cathode sensitivity with filters interposed (detailed below); the overall voltage for 200A/lm and the relevant dark current (at 20°C).  
A Corning glass filter (CS-5-58 ground to half stock thickness) is used to give a measure of the 'blue' sensitivity; a Corning glass filter (CS-2-62), which passes all radiation longer than 600nm, to indicate 'red' sensitivity, and a Wratten 87 filter, which passes all radiation longer than approximately 800nm, to indicate sensitivity in the near infra-red region.
- 2 Test data is obtained with K to d1 held at 150V and the 'standard' dynode chain.\*
- 3 Generally, tubes should be operated at, or near their rated overall sensitivity. Care should be taken not to exceed the maximum rated sensitivity or voltage.
- 4 For optimum stability under d.c. conditions, the mean anode current should not exceed 2 $\mu\text{A}$ .

\*For recommended dynode chains, refer to groups H, I, J on page 14 of the EMI Photomultiplier Tube Catalogue ref. P001/fP70 (available on request).

### Mechanical Characteristics

**Envelope diameter:** Maximum 29.00 mm (1.14 in)

**Cathode diameter:** Nominal 23.00 mm (0.91 in)

**Cathode type:** 9824 Bialkali (KCs) D125 Trialkali (S20)

**Window material:** 9824B, D125, Corning 9741 9824QB, D125QB quartz (spectrosil)

**Dynodes:** 11 box and grid with Cs Sb surfaces

**Base:** B14B low loss pressed glass (socket supplied)




### RATINGS

	9824	D125	
<b>Overall sensitivity: rated</b>	200A/lm	200A/lm	
	<b>maximum</b> 2000A/lm	2000A/lm	
<b>Voltage, cathode to d1: recommended maximum</b>	150V 300V	150V 300V	
<b>Voltage, anode to cathode: maximum</b>	2000V	2000V	
<b>Anode current (mean): maximum</b>	0.1mA	0.1mA	
<b>Anode dissipation: maximum</b>	0.1W	0.1W	
<b>Cathode current: maximum (using whole area)</b>	0.05 $\mu\text{A}$	1 $\mu\text{A}$	
<b>Anode pulse rise time: typical</b>	12 n. sec.	12 n. sec.	
<b>Anode pulse f.w.h.m.: typical</b>	50 n. sec.	50 n. sec.	
<b>Transit time: typical</b>	70 n. sec.	70 n. sec.	
<b>Capacitance, anode to all dynodes: typical</b>	6 pF	6 pF	
<b>Operating temperature: maximum</b>	60°C	60°C	
	<b>minimum</b>	-5°C	-180°C
<b>Dark current shot noise equivalent input*</b>	<b>Lumens</b>	$8.0 \times 10^{-14}$	$4.5 \times 10^{-13}$
	<b>Watts</b>	$6.5 \times 10^{-17}$	$8.4 \times 10^{-16}$


\* Calculated from typical performance data using Q.E. at  $\lambda$  peak and assuming  $\Delta f$  of 1 Hz and enhancement factor of unity.

Tube Type Number	Cathode Sensitivity				Overall Sensitivity						
	$\mu\text{A}/\text{lm}$		Corning Blue		V. Overall		2000A/lm		V. Overall		2000A/lm
	Min.	Typ.	Min.	Typ.	Typ.	Max.	Dark Current nA	Max.	Typ.	Max.	Dark Current nA
9824B 9824QB	—	50	5.0	7.5	1150	—	0.2	2.0	1550	—	2.0
D125B D125QB	80	120	—	—	1050	—	5.0	25	1450	—	50



SPECTRAL RESPONSE D125



SPECTRAL RESPONSE 9824



OVERALL SENSITIVITY vs OVERALL VOLTAGE



All dimensions are in millimetres with inches shown in parentheses.

# New 38mm (1½ in) diameter tubes having S20 or bialkali photocathode and 10 dynodes

## Provisional Data

- End-on low potassium borosilicate envelope with option of quartz (Spectrosil) window
- Fast response—3 n. sec. rise time
- Choice of plain glass (B14B) or overcapped (B12) base
- Bialkali photocathode (9633) for high blue sensitivity and very low dark current or trialkali (9683) for extended red

These new tubes are intended for use in applications where the high performance EMI 9635, 9750, 9558, 9658, 9659 tubes are not required. Typical applications include general low level photometry and laser detection.

Spectral coverage of the tubes with borosilicate window is about 310 to 650nm (bialkali) and about 310nm to beyond 850nm (trialkali). The U.V. response can be extended to approximately 170nm by using tubes with quartz windows.

### Notes

- 1 Each tube is individually calibrated and supplied with a test ticket giving the cathode sensitivity in  $\mu\text{A}/\text{lm}$  (except 9633); cathode sensitivity measurements with filters appropriate to the type of photocathode; the overall voltage for 20A/lm and the relevant dark current (at 20°C).  
A Corning glass filter (CS-5-58 ground to half stock thickness) is used to give a measure of the 'blue' sensitivity; a Corning glass filter (CS-2-62), which passes all radiation longer than approximately 600nm, to indicate 'red' sensitivity (except 9633), and a Wratten 87 filter, which passes all radiation longer than approximately 800 nm, to indicate sensitivity in the near infra red region (except 9633).
- 2 Test data is obtained with K to d1 held at 100 volts and the standard dynode chain.\*
- 3 Generally, tubes should be operated at, or near, their rated overall sensitivity. Care should be taken not to exceed the maximum rated sensitivity or voltage.
- 4 For optimum stability under d.c. conditions, the mean anode current should not exceed 1  $\mu\text{A}$ .

\*For recommended dynode chains, refer to groups D', E', F' on page 14 of the EMI Photomultiplier Tube Catalogue ref. P001/fP70 (available on request).

## MECHANICAL CHARACTERISTICS

Envelope diameter	Maximum	40 mm (1.58 in)
Cathode diameter	Nominal	25 mm (0.98 in)
Cathode type		9683 trialkali (S20) 9633 bialkali (KCs)
Window material		9633B, 9633KB, 9683B, 9683KB borosilicate






## RATINGS

	9633	9683
Overall sensitivity: rated	20A/lm	20A/lm
maximum	50A/lm	50A/lm
Voltage, cathode to d1: recommended	100V	100V
maximum	150V	150V
Voltage, anode to cathode: maximum	1600V	1600V
Anode current (mean): maximum	0.5mA	0.5mA
Anode dissipation: maximum	0.5W	0.5W
Cathode current: maximum (using whole area)	0.1 $\mu\text{A}$	1.5 $\mu\text{A}$
Anode pulse rise time: typical	2.5 n. sec.	2.5 n. sec.
Anode pulse f.w.h.m.: typical	6.0 n. sec.	6.0 n. sec.
Transit time: typical	25 n. sec.	25 n. sec.
Capacitance, anode to all dynodes: typical	7pF	7pF
Operating temperature: maximum	60°C	60°C
minimum	-5°C	-180°C
Dark current shot noise equivalent input*	Lumens $2.3 \times 10^{-13}$ Watts $2.3 \times 10^{-16}$	$4 \times 10^{-13}$ $1.4 \times 10^{-15}$

\* Calculated from typical performance data using Q.E. at  $\lambda$  peak and assuming  $\Delta f$  of 1 Hz and enhancement factor of unity.

	9633QB, 9633QKB, 9683QB, 9683QKB, quartz (Spectrosil)
Dynodes	10 compact focused with BeCu surfaces
Base	9633B, 9633QB, 9683B, 9683QB, B14B (socket supplied) 9633KB, 9633QKB, 9683KB, 9683QKB, B12 (socket extra)

Tube Type Number	Cathode Sensitivity				Overall Sensitivity							
	$\mu\text{A}/\text{lm}$		Corning Blue		20A/lm				50A/lm			
	Min.	Typ.	Min.	Typ.	V. Overall		Dark Current nA		V. Overall		Dark Current nA	
				Typ.	Max.	Typ.	Max.	Typ.	Max.	Typ.	Max.	
9633B	—	60	—	8.0	1300	—	0.2	—	1600	—	0.5	
9633QB	—	200	—	—	1250	—	2.0	—	1500	—	5.0	



All dimensions are in millimetres with inches shown in parentheses.

**PIN CONNECTIONS** (viewed from below starting left of short pin or key)

Type	1	2	3	4	5	6	7	8	9	10	11	12	13	14	socket
9633	d1	d3	d5	d7	d9	A	d10	d8	d6	d4	d2	—	K	—	B14A
9683															
9633K	d1	d3	d5	d7	d9	A	d10	d8	d6	d4	d2	K	—	—	B12A
9683K															

# New 50mm (2in) diameter tubes with trialkali (S20) photocathode and 6 dynodes

## Provisional Data

- Wide spectral coverage—310nm to 850nm (170 to 850nm quartz version)
- Small size—113 mm overall length
- Low dark current
- Surface of triangular prisms on inside of end window to enhance quantum efficiency.

This tube is a variant of the EMI 9658 type having a shorter envelope and 6 dynodes. The main application will be in the detection of relatively high light levels such as from laser or cathode ray tube sources, or generally where low gain is beneficial. The small number of dynodes also reduces the rate of change in gain with voltage.

An overcapped version is available, which is designated D214KB.

### Notes

- 1 Each tube is individually calibrated and supplied with a test ticket giving the cathode sensitivity in  $\mu\text{A}/\text{lm}$ ; cathode sensitivity with filters interposed (detailed below); the overall voltage for 1A/lm and the relevant dark current (at 20°C).  
A Corning glass filter (CS-5-58 ground to half stock thickness) is used to give a measure of the 'blue' sensitivity; a Corning glass filter (CS-2-62), which passes all radiation longer than 600 nm, to indicate 'red' sensitivity, and a Wratten 87 filter, which passes all radiation longer than approximately 800nm, to indicate sensitivity in the near infra-red region.
- 2 Test data is obtained with K to d1 held at 150V and a 'linear' dynode chain.\*
- 3 Generally, tubes should be operated at, or near their rated overall sensitivity. Care should be taken not to exceed the maximum rated sensitivity or voltage.
- 4 For optimum stability under d.c. conditions, the mean anode current should not exceed 10 $\mu\text{A}$ .

\*For recommended dynode chains, refer to groups A, B, C on page 14 of the EMI Photomultiplier Tube Catalogue ref. P001/fP70 (available on request).

### MECHANICAL CHARACTERISTICS

Envelope diameter	Maximum	51.5 mm (2.02 in)
Cathode diameter	Nominal	40.0 mm (1.57 in)
Cathode type	Trialkali (S20)	
Window material	Borosilicate	
Dynodes	6 venetian blind with CsSb surfaces	
Base	B19A low loss pressed glass (socket supplied)	




### RATINGS

Overall sensitivity: rated maximum	1A/lm 5A/lm	
Voltage, cathode to d1: recommended maximum	150V 300V	
Voltage, anode to cathode: maximum	2000V	
Anode current (mean): maximum	1mA	
Anode dissipation: maximum	1W	
Cathode current: maximum (using whole area)	5 $\mu\text{A}$	
Anode pulse rise time: typical	10 n. sec.	
Anode pulse f.w.h.m.: typical	14 n. sec.	
Transit time: typical	40 n. sec.	
Capacitance, anode to all dynodes: typical	8pF	
Operating temperature: maximum minimum	60°C -180°C	
Dark current shot noise equivalent input*	Lumens 'Watts	0.8 × 10 <sup>-13</sup> 4.3 × 10 <sup>-16</sup>



\* Calculated from typical performance data using Q.E. at  $\lambda$  peak and assuming  $\Delta f$  of 1 Hz and enhancement factor of unity.

Tube Type Number	Cathode Sensitivity $\mu\text{A}/\text{lm}$	Overall Sensitivity					
		1A/lm			5A/lm		
		V. Overall	Dark Current nA	V. Overall	Dark Current nA	V. Overall	Dark Current nA
D214B	Min. Typ.	Typ. Max.	Typ. Max.	Typ. Max.	Typ. Max.	Typ. Max.	Typ. Max.
	200 350	800 1800	0.5 5.0	1100 —	2.5		

OVERALL SENSITIVITY vs OVERALL VOLTAGE



SPECTRAL RESPONSE



All dimensions in millimetres (inches in parenthesis).

**PIN CONNECTIONS** (viewed from below starting left of short pin or key)

Type	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	socket	
D214B	d1	—	—	—	—	—	—	A	—	d6	d5	d4	d3	d2	—	—	—	—	—	K	B19A
D214KB	i.c.	i.c.	i.c.	i.c.	i.c.	A	i.c.	d6	d5	d4	d3	d2	K	d1	—	—	—	—	—	—	B14A

i.c.—internal connection

# New 50mm (2in) diameter tube with trialkali(S20) photocathode and 11 dynodes

- High sensitivity in the blue and green parts of the spectrum
- S20 photocathode which has high conductivity and freedom from "lag" effects
- Physically interchangeable with other 50 mm tube in EMI S20 range


The tube is an addition to the existing range of EMI 50 mm tubes with S20 photocathodes. It is intended for use in the blue or green channels of colour flying spot equipment where high signal to noise ratio, under high light level conditions, is important. The main application is in telecine machines for broadcast use.

Depending upon availability, the 9658F may be supplied as an alternative. This tube meets the same performance specification and has identical external dimensions and pin connections.

## Notes

- 1 Each tube is individually calibrated and supplied with a test ticket giving the cathode sensitivity in  $\mu\text{A}/\text{lm}$ ; cathode sensitivity with filters interposed (detailed below); the overall voltage for 200A/lm and the relevant dark current (at 20°C).  
A Corning glass filter (CS-5-58 ground to half stock thickness) is used to give a measure of the 'blue' sensitivity; a Corning glass filter (CS-2-62), which passes all radiation longer than 600nm, to indicate 'red' sensitivity, and a Wratten 87 filter, which passes all radiation longer than approximately 800nm, to indicate sensitivity in the near infra-red region. Suitability for use in the blue or green channel depends upon a comparison between Corning red and Corning blue filter readings and tube tickets are marked accordingly.
- 2 Test data is obtained with K to d1 held at 150V and the 'standard' dynode chain.\*
- 3 Generally, tubes should be operated at, or near, their rated overall sensitivity. Care should be taken not to exceed the maximum rated sensitivity or voltage.

\*For recommended dynode chains, refer to groups H, I, J on page 14 of the EMI Photomultiplier Tube Catalogue ref. P001/fP70 (available on request).



## RATINGS




<b>Overall sensitivity: rated maximum</b>	<b>200A/lm</b>
	<b>2000A/lm</b>
<b>Voltage, cathode to d1: recommended maximum</b>	<b>300V</b>
	<b>350V</b>
<b>Voltage, anode to cathode: maximum</b>	<b>2200V</b>
<b>Anode current (mean): maximum</b>	<b>1mA</b>
<b>Anode dissipation: maximum</b>	<b>1W</b>
<b>Cathode current: maximum (using whole area)</b>	<b>5<math>\mu</math>A</b>
<b>Anode pulse rise time: typical</b>	<b>10 n. sec.</b>
<b>Anode pulse f.w.h.m.: typical</b>	<b>15 n. sec.</b>
<b>Transit time: typical</b>	<b>55 n. sec.</b>
<b>Capacitance, anode to all dynodes: typical</b>	<b>8pF</b>
<b>Operating temperature: maximum minimum</b>	<b>60°C</b>
	<b>-180°C</b>
<b>Dark current shot noise equivalent input*</b>	<b>Lumens 1.5 <math>\times 10^{-13}</math></b>
	<b>Watts 3.7 <math>\times 10^{-16}</math></b>

\* Calculated from typical performance data using Q.E. at  $\lambda$  peak and assuming  $\Delta f$  of 1 Hz and enhancement factor of unity.


## MECHANICAL CHARACTERISTICS

<b>Envelope diameter</b>	<b>Maximum 51.5 mm (2.02 in)</b>
<b>Cathode diameter</b>	<b>Nominal 9558F 44 mm (1.73 in)</b> <b>(9658F 40 mm)</b>
<b>Cathode type</b>	<b>Trialkali (S20)</b>
<b>Window material</b>	<b>Borosilicate</b>
<b>Dynodes</b>	<b>11 Venetian blind with CsSb surfaces</b>
<b>Base</b>	<b>B19A low loss pressed glass (socket supplied)</b>

Tube Type Number	Cathode Sensitivity		Overall Sensitivity								
			200A/lm				2000A/lm				
	Min.	Typ.	Corning Blue		V. Overall		Dark Current nA		V. Overall		Dark Current nA
		Min.	Typ.	Typ.	Max.	Typ.	Max.	Typ.	Max.	Typ.	
9558F	100	150	7.0	9.0	1100	2000	5.0	50	1450	—	50



All dimensions are in millimetres with inches shown in parentheses



**Pin Connections Types 9558, 9658 B19A Socket (Viewed from below)**

# New 50mm (2in) diameter tubes with bialkali photocathode and 11 or 13 dynodes

- High 'blue' response and very low dark current
- 9804B high gain, typically  $4 \times 10^7$  at 1000V
- 9757B medium gain, typically  $4 \times 10^6$  at 950V

The 9757B and 9804B are intended to replace the 6097S and 9514S types. Being physically interchangeable, they offer the advantage of an immediate improvement in blue response together with a dark current level as low as that obtained from the 'S' type tubes. The venetian blind dynode structures utilise highly efficient and stable CsSb secondary emitting surfaces.

The main applications for the new tubes are in low level photometry, nuclear radiation monitoring and thermoluminescent dosimetry.

## Notes

- 1 Each tube is individually calibrated and supplied with a test ticket giving the cathode sensitivity in  $\mu\text{A}/\text{lm}$ ; cathode sensitivity measured with a Corning blue CS-5-58 filter (polished to half stock thickness) interposed; the overall voltage for 200A/lm (9757) or 2000A/lm (9804) and the relevant dark current (at 20 °C). The quantum efficiency at 420nm can be estimated by multiplying the Corning blue reading by 2.5.
- 2 Test data is obtained with K to d1 held at 150 volts and the 'standard' dynode chain.\*
- 3 Generally, tubes should be operated at, or near, their rated overall sensitivity. Care should be taken not to exceed the maximum rated sensitivity or voltage.
- 4 For optimum stability under d.c. conditions, the mean anode current should not exceed 10 $\mu\text{A}$ .

\*For recommended dynode chains, refer to groups H, I, J, (9757) or K, L, M (9804) on page 14 of the EMI Photomultiplier Catalogue ref. P001/fP70 (available on request).

## MECHANICAL CHARACTERISTICS

Envelope diameter	Maximum	51.5 mm (2.02 in)
Cathode diameter	Nominal	44 mm (1.73 in)
Cathode type	Bialkali (SbKCs)	
Window material	Borosilicate	
Dynodes	9757 (11 stages); 9804 (13 stages); venetian blind dynodes with CsSb secondary emitting surfaces	
Base	Low loss 15-pin pressed glass base furnished with high quality Teflon socket type B15B	



## RATINGS

	9757	9804	
Overall sensitivity: rated maximum	200A/lm 2000A/lm	2000A/lm 5000A/lm	
Voltage, cathode to d1: recommended maximum	150V 300V	150V 300V	
Voltage, anode to cathode: maximum	1800V	2500V	
Anode current (mean): maximum	0.1mA	0.1mA	
Anode dissipation: maximum	0.1W	0.1W	
Cathode current: maximum (using whole area)	0.1 $\mu\text{A}$	0.1 $\mu\text{A}$	
Anode pulse rise time: typical	11 n. sec.	11 n. sec.	
Anode pulse f.w.h.m.: typical	15 n. sec.	16 n. sec.	
Transit time: typical	55 n. sec.	60 n. sec.	
Capacitance, anode to all dynodes: typical	8pF	8pF	
Operating temperature: maximum minimum	60°C -5°C	60°C -5°C	
Dark current shot noise equivalent input*	Lumens Watts	$8.0 \times 10^{-14}$ $4.8 \times 10^{-17}$	$8.0 \times 10^{-14}$ $4.8 \times 10^{-17}$


\* Calculated from typical performance data using Q.E. at  $\lambda$  peak and assuming  $\Delta f$  of 1 Hz and enhancement factor of unity.

Tube Type Number	Cathode Sensitivity		Overall Sensitivity								
			200A/lm				2000A/lm				
	$\mu\text{A}/\text{lm}$	Corning Blue	V. Overall		Dark Current nA		V. Overall		Dark Current nA		
Min.	Typ.	Min.	Typ.	Typ.	Max.	Typ.	Max.	Typ.	Max.	Typ.	
9757B	—	50	5.0	9.0	950	1400	0.2	1.0	1250	—	2.0
9804B	—	50	5.0	9.0	1000	1700	2.0	10	1100	—	5.0


The maximum dark current count rate is 1000 electrons/second.

OVERALL SENSITIVITY vs OVERALL VOLTAGE

9757





SPECTRAL RESPONSE




OVERALL SENSITIVITY vs OVERALL VOLTAGE

9804



All dimensions are in millimetres with inches shown in parentheses



Pin Connections Types 9757\* 9804

B15B Socket (Viewed from below)

\*Pins 5&7 IC

# New 13 stage 50mm (2in) diameter tubes with bialkali photocathode having an effective diameter of 10mm

- High 'blue' response with extremely low dark current
- High gain, typically  $4 \times 10^7$  at 1150V
- Choice of borosilicate (310nm to 650nm) or quartz (170nm to 650nm) windows

This tube is intended to replace the 9502S and 6256S types with which it is physically interchangeable. It has a unique K-d<sub>1</sub> geometry, together with a venetian blind dynode structure having highly stable and efficient CsSb secondary emitting surfaces. This combination of dynode system and bialkali photocathode makes the 9789 capable of giving very high values of gain together with extremely low dark currents.

Main applications include very low level photometry, photon counting and astronomy.

## Notes

- 1 Each tube is individually calibrated and supplied with a test ticket giving the cathode sensitivity in  $\mu\text{A}/\text{lm}$ ; cathode sensitivity measured with a Corning blue CS-5-58 filter (polished to half stock thickness) interposed; the overall voltage for 200A/lm and the relevant dark current (at 20°C). The quantum efficiency at 420nm can be estimated by multiplying the Corning blue reading by 2.5.
- 2 Test data is obtained with K to d1 held at 150 volts and the 'standard' dynode chain.\*
- 3 Generally, tubes should be operated at, or near, their rated overall sensitivity. Care should be taken not to exceed the maximum rated sensitivity or voltage.
- 4 For optimum stability under d.c. conditions, the mean anode current should not exceed  $10\mu\text{A}$ .

\*For recommended dynode chains, refer to groups K, L, M, on page 14 of the EMI Photomultiplier Catalogue ref. P001/fP70 (available on request).



## RATINGS

<b>Overall sensitivity: rated maximum</b>	<b>2000A/lm 5000A/lm</b>
<b>Voltage, cathode to d1: recommended maximum</b>	<b>150V 300V</b>
<b>Voltage, anode to cathode: maximum</b>	<b>2100V</b>
<b>Anode current (mean): maximum</b>	<b>1mA</b>
<b>Anode dissipation: maximum</b>	<b>1W</b>
<b>Cathode current: maximum (using whole area)</b>	<b>0.3<math>\mu\text{A}</math></b>
<b>Anode pulse rise time: typical</b>	<b>10 n. sec.</b>
<b>Anode pulse f.w.h.m.: typical</b>	<b>15 n. sec.</b>
<b>Transit time: typical</b>	<b>55 n. sec.</b>
<b>Capacitance, anode to all dynodes: typical</b>	<b>8pF</b>
<b>Operating temperature: maximum minimum</b>	<b>60°C -5°C</b>
<b>Dark current shot noise equivalent input*</b>	<b>Lumens Watts 2.5 <math>\times 10^{-14}</math> 1.5 <math>\times 10^{-17}</math></b>

\* Calculated from typical performance data using Q.E. at  $\lambda$  peak and assuming  $\Delta f$  of 1 Hz and enhancement factor of unity.


## MECHANICAL CHARACTERISTICS

<b>Envelope diameter</b>	<b>Maximum 51.5 mm (2.02 in)</b>	<b>Dynodes</b>	<b>13 venetian blind dynodes with CsSb secondary emitting surfaces</b>
<b>Cathode diameter</b>	<b>Minimum 10 mm (0.39 in)</b>		
<b>Cathode type</b>	<b>Bialkali (SbKCs)</b>		
<b>Window material</b>	<b>9789B Borosilicate, 9789QB Spectrosil (fused silica)</b>	<b>Base</b>	<b>Low loss 15-pin pressed glass furnished with high quality Teflon socket type B15B</b>




Tube Type Number	Cathode Sensitivity		Overall Sensitivity								
	$\mu\text{A}/\text{lm}$		2000A/lm				5000A/lm				
	Min.	Typ.	Corning Blue		V. Overall		Dark Current nA		V. Overall		Dark Current nA
			Min.	Typ.	Typ.	Max.	Typ.	Max.	Typ.	Max.	Typ.
9789B	—	50	5.0	8.0	1150	1700	0.2	1.0	1300	—	0.5
9789QB	—	50	5.0	8.0	1150	1700	0.2	1.0	1300	—	0.5

The maximum dark current count rate is 50 electrons/second.

OVERALL SENSITIVITY vs OVERALL VOLTAGE



SPECTRAL RESPONSE



Pin Connections Type

9789

B15B Socket

(Viewed from below)

All dimensions are in millimetres with inches shown in parentheses

# New 75mm (3in) diameter tubes with bialkali photocathode and an internal end window surface of triangular prisms

## Provisional Data

- Very high 'blue' response and low dark current
- Good energy resolution, e.g. 12.5% f.w.h.m. to  $^{57}\text{Co}$
- Surface of triangular prisms on inside of end window

The 9778 is a new tube based on the EMI 9758 with which it is physically interchangeable. However the 9778 has an inside end window surface which promotes multiple internal reflection of the incoming light. This gives a better utilisation of the incident light but its main purpose for use in gamma cameras is to eliminate abrupt changes of sensitivity with varying angles of incidence, which exist in normal semi-transparent tubes with plane windows. The tube is suitable for use in gamma cameras and other low level scintillation counters, gamma and X-ray spectrometers. Tubes can be supplied in batches having a narrow spread of voltages for 50A/lm overall sensitivity.

### Notes

- 1 Each tube is individually calibrated and supplied with a test ticket giving the cathode sensitivity in  $\mu\text{A}/\text{lm}$ ; cathode sensitivity measured with a Corning blue CS-5-58 filter (polished to half stock thickness) interposed; the overall voltage for 50A/lm and the relevant dark current (at 20°C). The quantum efficiency at 420nm can be estimated by multiplying the Corning blue reading by 2.5.
- 2 Test data is obtained with K to d1 held at 300 volts and the 'standard' dynode chain.\*
- 3 Generally, tubes should be operated at, or near, their rated overall sensitivity. Care should be taken not to exceed the maximum rated sensitivity or voltage.
- 4 For optimum stability under d.c. conditions, the mean anode current should not exceed 10 $\mu\text{A}$ .

\*For recommended dynode chains, refer to groups D, E, F on page 14 of the EMI Photomultiplier Catalogue ref. P001/fP70 (available on request).

## MECHANICAL CHARACTERISTICS

Envelope diameter	Maximum	78 mm (3.07 in)
Cathode diameter	Nominal	65 mm (2.56 in)
Cathode type		Bialkali (KCs)
Window material		Borosilicate (internal surface consists of 3 mm square base triangular prisms)



## RATINGS


Overall sensitivity: rated maximum	50A/lm 200A/lm	
Voltage, cathode to d1: recommended maximum	300V 350V	
Voltage, anode to cathode: maximum	2100V	
Anode current (mean): maximum	1mA	
Anode dissipation: maximum	1W	
Cathode current: maximum (using whole area)	0.5 $\mu\text{A}$	
Anode pulse rise time: typical	12 n. sec.	
Anode pulse f.w.h.m.: typical	50 n. sec.	
Transit time: typical	70 n. sec.	
Capacitance, anode to all dynodes: typical	8pF	
Operating temperature: maximum minimum	60°C -5°C	
Dark current shot noise equivalent input*	Lumens Watts	2 $\times 10^{-13}$ 2 $\times 10^{-16}$

\* Calculated from typical performance data using Q.E. at  $\lambda$  peak and assuming  $\Delta f$  of 1 Hz and enhancement factor of unity.


Dynodes	9 venetian blind with CsSb surfaces
Base	9778B B19A low loss pressed glass (socket supplied) 9778KB overcapped B14A (socket extra)

Tube Type Number	Cathode Sensitivity		Overall Sensitivity								
			50A/lm				200A/lm				
	$\mu\text{A}/\text{lm}$		Corning Blue		V. Overall		Dark Current nA		V. Overall		Dark Current nA
Min.	Typ.	Min.	Typ.	Typ.	Max.	Typ.	Max.	Typ.	Max.	Typ.	
9778B		85	7.0	11.0	1050	1750	0.5	50	1250	—	2.0



OVERALL SENSITIVITY vs OVERALL VOLTAGE



SPECTRAL RESPONSE



RELATIVE SENSITIVITY VARIATION WITH CHANGING ANGLE OF INCIDENCE



All dimensions are in millimetres with inches shown in parentheses.

**PIN CONNECTIONS** (viewed from below starting left of short pin or key)

Type	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	socket
9778B	i.c.	d1	d2	d3	d4	d5	d6	d7	d8	i.c.	d9	A	i.c.	i.c.	i.c.	i.c.	i.c.	i.c.	K	B19A
9778KB	d1	d2	d3	d4	d5	d6	d7	d8	i.c.	d9	A	i.c.	i.c.	K	—	—	—	—	—	B14A

i.c.—internal connection

# New 50 mm (2in) diameter tube with bialkali photocathode on a curved end window

## Provisional Data

- High blue sensitivity and low dark current
- Curved end window with matt surface on outer face for coupling to liquid scintillation spectrometer reflective chambers
- Low background from low potassium glass

The 9820 is a new tube based on the EMI 9813, which has a linear focused dynode system for fast response and a high efficiency bialkali photocathode. The tube is intended for liquid scintillation counting in coincidence arrangements and will give a high  $^3\text{H}$  efficiency with a low background. The use of low potassium glass for the envelope and window ensures that the background associated with the tube is kept to a minimum. Tubes can be supplied in pairs for use in coincidence systems.

### Notes

- 1 Each tube is individually calibrated and supplied with a test ticket giving the cathode sensitivity in  $\mu\text{A}/\text{lm}$ ; cathode 'blue' sensitivity measured with a Corning glass filter (CS-5-58 ground to half stock thickness) from which the quantum efficiency at 420nm can be estimated by multiplying by 2.50; the overall voltage for 200A/lm and the relevant dark current (at 20°C).
- 2 Test data is obtained with K to d1 held at 300 volts and a linear test dynode chain.
- 3 F (focus) should be connected to d1 for normal operation.
- 4 Generally, tubes should be operated at or near their rated overall sensitivity. Care should be taken not to exceed the maximum rated sensitivity or voltage.
- 5 For optimum stability under d.c. conditions, the mean anode current should not exceed 10 $\mu\text{A}$ .
- 6 The glass envelope is coated with a screen which is connected to the cathode pin and over-wound with 3M's electrical tape, to ensure optimum dark current performance. Special care should be taken if the tube is operated with the cathode at high negative potential.

### DYNODE CHAIN DESIGN

For best results a linear dynode chain should be used. The circuit shown on page 26 is recommended but with resistance R used between each pair of stages from d3 onwards.



### RATINGS

<b>Overall sensitivity: rated maximum</b>	<b>200A/lm</b>	<b>2000A/lm</b>
<b>Voltage, cathode to d1: recommended maximum</b>	<b>300V</b>	<b>500V</b>
<b>Voltage between dynodes: maximum</b>	<b>550V</b>	
<b>Voltage, anode to last dynode: maximum</b>	<b>550V</b>	
<b>Voltage, anode to cathode: maximum</b>	<b>2500V</b>	
<b>Anode current (mean) maximum</b>	<b>0.1mA</b>	
<b>Anode dissipation: maximum</b>	<b>0.1W</b>	
<b>Cathode current: maximum (using whole area)</b>	<b>0.3<math>\mu\text{A}</math></b>	
<b>Anode pulse rise time: typical</b>	<b>2.0 n. sec.</b>	
<b>Anode pulse f.w.h.m.: typical</b>	<b>2.9 n. sec.</b>	
<b>Transit time: typical</b>	<b>40 n. sec.</b>	
<b>Capacitance, anode to all dynodes: typical</b>	<b>5pF</b>	
<b>Operating temperature: maximum minimum</b>	<b>60°C</b>	<b>-20°C</b>
<b>Dark current shot noise equivalent input*</b>	<b>Lumens</b>	<b>1.6 <math>\times 10^{-13}</math></b>
	<b>Watts</b>	<b>1.3 <math>\times 10^{-16}</math></b>




\* Calculated from typical performance data using Q.E. at  $\lambda$  peak and assuming  $\Delta f$  of 1 Hz and enhancement factor of unity.

### MECHANICAL CHARACTERISTICS


<b>Envelope diameter</b>	<b>Maximum 53 mm (2.09 in)</b>
<b>Cathode diameter</b>	<b>Nominal 46 mm (1.81 in)</b>
<b>Cathode type</b>	<b>Bialkali (KCsSb)</b>
<b>Window material</b>	<b>Borosilicate or Spectrosil</b>
<b>Dynodes</b>	<b>12 linear focused with CsSb surfaces</b>
<b>Base</b>	<b>B19A low loss pressed glass (socket supplied)</b>

Tube Type Number	Cathode Sensitivity				Overall Sensitivity $\phi$							
	$\mu\text{A/lm}$		Corning Blue		200A/lm				2000A/lm			
	Min.	Typ.	Min.	Typ.	V. Overall		Dark Current nA		V. Overall		Dark Current nA	
				Typ.	Max.	Typ.	Max.	Typ.	Max.	Typ.	Max.	Typ.
9820B QB	—	65	7.0	10.0	1400	2000	1.0	4.0	1800	—	—	10

$\phi$  measured using linear dynode chain.



All dimensions in millimetres (inches in parenthesis).



# 50 mm (2in) diameter fast linear focused tubes with S11, bialkali, or S20 photocathodes

- Semi-transparent curved photocathode surface on end window which is flat externally
- Linear focused dynode structure with caesiated Be O surfaces for high gain and good stability
- Spectral responses S-11, S-13, Bialkali or S-20
- 10, 12 or 14 dynode stages
- Flexibility of envelope configuration, i.e. B19A(uncapped) or B20-102 (overcapped) base; borosilicate (Pyrex) or quartz (Spectrosil) end window material

The 2 in. diameter tubes 9810-9818 are designed for applications requiring very fast response times such as the detection and measurement of short lived nuclear events and pulsed light sources.

This new series has been introduced to replace the 9594 range. Apart from being more comprehensive, it incorporates design modifications to give improved reliability and performance. The external dimensions and pin connections are unchanged to preserve interchangeability.



## MECHANICAL CHARACTERISTICS

Maximum envelope diameter	All types 53 mm (2.09 in)
Nominal cathode diameter	All types 46 mm (1.81 in)
Window material (see drawings for shapes)	'B' and 'KB' types borosilicate (Pyrex); 'QB' and 'QKB' types quartz (Spectrosil)
Dynodes	Caesiated Be O surface, linear focused
Base	'B' and 'QB' types B19A (socket supplied) 'KB' and 'QKB' types Jedec B20-102 (socket extra)


## PRELIMINARY SELECTION GUIDE

Basic Type*	Spectral Response	Number of Dynodes	Old Type	Remarks
9810B	S-11	14	9594UB	General purpose high gain parent type
9811B	S-11	12	—	Medium gain version of 9810B
9812B	S-11	10	9595UB	Lower gain version of 9810B for higher light levels
9813B	Bialkali	14	—	High gain, low dark current, high "blue" response
9814B	Bialkali	12	—	Medium gain version of 9813B
9815B	Bialkali	10	—	Lower gain version of 9813B for higher light levels
9816B	S-20	14	9597UB	High gain, spectral response extends to IR
9817B	S-20	12	—	Medium gain version of 9816B
9818B	S-20	10	9596UB	Lower gain version of 9816B, for higher light levels

\*Versions with quartz end window and/or overcapped are available in which case Q and/or K respectively are added to the type number, e.g. 9813QKB.

**Notes**

- 1 a) Each tube is individually calibrated and supplied with a test ticket giving the cathode sensitivity in  $\mu\text{A}/\text{lm}$ ; cathode sensitivity measurements with filters appropriate to type of photocathode; the overall voltages for 50A/lm (10-stage tubes only), 200 A/lm (12-stage tubes only), 5000 A/lm (14-stage tubes only), and the relevant dark current (at 20°C).  
A Corning glass filter (CS-5-58 ground to half stock thickness) is used to give a measure of the "blue" sensitivity; a Corning glass filter (CS-2-62), which passes all radiation longer than approximately 600nm, to indicate "red" sensitivity, and a Wratten 87 filter, which passes all radiation longer than approximately 800nm, to indicate sensitivity in the near infra-red region.
- b) Test data is obtained with K to d1 held at 300 volts and the non-linear test dynode chain.
- c) F (focus) should be connected to d1 for normal operation.
- d) Generally, tubes should be operated at or near their rated overall sensitivity. Care should be taken not to exceed the maximum rated sensitivity or voltage.
- 2 For optimum stability under d.c. conditions, the mean anode current should not exceed  $10\mu\text{A}$ .
- 3 The glass envelope is coated with a screen which is connected to the cathode pin and over-wound with 3M's electrical tape, to ensure optimum dark current performance. Special care should be taken if the tube is operated with the cathode at high negative potential.
- 4 Overall sensitivity is greatly affected by variations in overall voltage, particularly in the case of tubes with a large number of dynodes.
- 5 In order to obtain a high value of peak output current, it is necessary to provide high voltages between the last dynode stages. The EMI test dynode chain is recommended for optimum results.




Tube type No.*	Cathode Sensitivity $\mu\text{A}/\text{lm}$		Overall Sensitivity 5000A/lm				Overall Sensitivity 10,000A/lm	
	Min.	Typ.	V Overall Typ.	Dark Current nA Typ.	Max.	V Overall Typ.	Dark Current nA Typ.	
9810B†	50	75	2250	2700	100	1000	2380	200
9810B‡			1980	—			2080	
9813B†	§	§	2350	2700	10	100	2480	20
9813B‡			2050	—			2160	
9816B†	100	170	2150	2700	50	1000	2270	100
9816B‡			1900	—			2000	
			Overall Sensitivity 200A/lm				Overall Sensitivity 2000A/lm	
9811B†	50	75	2050	2500	15	100	2600	150
9811B‡			1750	—			2220	
9814B†	§	§	2150	2500	2	8	2750	20
9814B‡			1850	—			2350	
9817B†	100	170	1950	2500	10	100	2470	100
9817B‡			1650	—			2070	
			Overall Sensitivity 50A/lm				Overall Sensitivity 500A/lm	
9812B†	50	75	1850	2300	2	50	2500	20
9812B‡			1560	—			2100	
9815B†	§	§	1950	2300	1	5	2620	10
9815B‡			1640	—			2200	
9818B†	100	170	1750	2300	2	50	2350	20
9818B‡			1470	—			1970	

\* Tube characteristics apply equally to overcapped and quartz versions.  
† Non-linear dynode chain.  
‡ Linear dynode chain.  
§ Types 9813, 9814 9815 have a minimum Corning blue cathode sensitivity of 7.0 and typical of 10.0.

**DYNODE CHAIN DESIGN**

The EMI (Non-linear) test dynode chain is shown below :



$R' = 1.4R$  for normal operation but may be varied between  $R$  and  $1.4R$  without significant change in performance  
 $RL = 50 \Omega$   $n =$  number of dynodes  
 $Z = \text{IS 4150A (150V)}$   $c = 1000\text{pF}$

The value of  $R$  should be chosen to provide a series current, in the divider chain, of at least 10 times the mean expected anode current. To maintain linearity when high current anode pulses are being drawn, the voltages supplying the final dynodes should be increased, and if resistors are being used they should be decoupled by suitably large value capacitors. The photomultiplier output socket will also affect operation under pulsed conditions.

General notes on design aspects are given in the Introduction of the EMI photomultiplier catalogue P001/fP70, available on request.

## RATINGS

PARAMETER	TUBE TYPE (All versions)		
	9810 9813 9816	9811 9814 9817	9812 9815 9818
Overall sensitivity: Rated	5,000 A/lm	200 A/lm	50 A/lm
Maximum	10,000 A/lm	2,000 A/lm	500 A/lm
Anode pulse rise time: Typical	2.4	2.0	1.6
(n. sec) Maximum	3.0	2.5	2.2
Anode Pulse f.w.h.m.: Typical	3.6	2.8	2.5
(n. sec) Maximum	4.4	4.0	3.5
Electron transit time: Typical	45	39	32
(n. sec)			
Voltage cathode to d1: Recommended	300	300	300
Maximum	500	500	500
Voltage between dynodes: Maximum	550	550	550
Voltage anode to last dynode: Maximum	550	550	550
Voltage anode to cathode: Maximum	3000	2800	2500
Maximum anode current (mean):	0.2mA	0.2mA	0.1mA
Maximum anode dissipation:	0.1W	0.1W	0.05W
Maximum tolerable cathode current (assuming whole cathode used)	S-11 types S-20 types Bialkali types		0.3 μA 5 μA 0.3 μA
Maximum operating temperature:	60°C	60°C	60°C
Minimum operating temperature:	-80°C*	-80°C*	-80°C*
Capacitance to all electrodes: Anode	B types 5pf	KB types 7pf	
Last dynode	B types 7pf	KB types 9pf	
Energy resolution using Na I (TI) crystal measured on type 9813	f.w.h.m. typically 7.5% for <sup>137</sup> Cs P/V typically 6.0:1 for <sup>60</sup> Co		

\* -20°C for 9813, 9814 and 9815

## DARK CURRENT SHOT NOISE EQUIVALENT INPUT\*

Type	9810	9811	9812	9813	9814	9815	9816	9817	9818
Lumens	$2.9 \times 10^{-13}$	$5.7 \times 10^{-13}$	$4.1 \times 10^{-13}$	$9.5 \times 10^{-14}$	$2.1 \times 10^{-13}$	$3.0 \times 10^{-13}$	$1.4 \times 10^{-13}$	$3.1 \times 10^{-13}$	$2.7 \times 10^{-13}$
Watts	$3.6 \times 10^{-16}$	$6.9 \times 10^{-16}$	$5.0 \times 10^{-16}$	$8.1 \times 10^{-17}$	$1.9 \times 10^{-16}$	$2.5 \times 10^{-16}$	$3.4 \times 10^{-16}$	$7.7 \times 10^{-16}$	$2.0 \times 10^{-15}$

\*Calculated from typical performance data using Q.E. at  $\lambda$  peak and assuming  $\Delta f$  of 1Hz and enhancement factor of unity.

## TUBE REPLACEMENT GUIDE

RCA	EMI*	PHILIPS	EMI†
4459	9817KB	56AVP	9810KB
6810	9810KB	56TUVF	9816QKB
6810A	9810KB	56DVP	9813KB
7264	9810KB	56DUVP	9813QKB
7265	9816KB	PHILIPS	EMI*
7326	9818KB	XP1000	9812KB
7746	9812KB	XP1001	9812KB
7850	9811KB	XP1002	9818KB
8575	9814QB	XP1003	9818QKB

## NOTES ON TIME CHARACTERISTICS

The anode pulse rise time is the time taken for the output pulse to rise from 10% to 90% of the peak when the photocathode of the tube is illuminated by a flash of light of very short duration.

The anode pulse f.w.h.m. is the full width of the output pulse measured at half maximum amplitude.

The electron transit time is the time difference between the arrival of a flash of light at the photocathode and the instant when the output pulse is a maximum.


In each case the whole photocathode is illuminated by a point source situated about 50 mm above its centre. Measurements referred to in the "electrical ratings" were made at the maximum rated voltage for 10 and 12 stage tubes and the standard rated voltage for 14 stage tubes.

\*nearest equivalent given - see data for differences

†direct plug-in replacement (focus voltage will need adjustment)


**OVERALL SENSITIVITY vs OVERALL VOLTAGE**

Linear dynode chain 14 stage tubes 9810 9813 9816




**OVERALL SENSITIVITY vs OVERALL VOLTAGE**

Non-linear dynode chain 14-stage tubes 9810 9813 9816




**OVERALL SENSITIVITY vs OVERALL VOLTAGE**

Linear dynode chain 12 stage tubes 9811 9814 9817




**OVERALL SENSITIVITY vs OVERALL VOLTAGE**

Non-linear dynode chain 12 stage tubes 9811 9814 9817






### OVERALL SENSITIVITY vs OVERALL VOLTAGE

Linear dynode chain 10 stage tubes 9812 9815 9818



### OVERALL SENSITIVITY vs OVERALL VOLTAGE

Non-linear dynode chain 10 stage tubes 9812 9815 9818



All dimensions in millimetres (inches in parenthesis).

#### PIN CONNECTIONS (viewed from below starting left of short pin or key)

Type	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	SOCKET
<b>9810B 9813B 9816B</b>	ic	d1	d3	d5	d7	d9	d11	d13	A	d14	d12	d10	d8	d6	d4	d2	ic	F	K	—	<b>B19A</b>
<b>9811B 9814B 9817B</b>	ic	d1	d1	d3	d5	d7	d9	d11	A	d12	d10	d8	d6	d4	d2	ic	ic	F	K	—	<b>B19A</b>
<b>9812B 9815B 9818B</b>	ic	d1	d1	d1	d3	d5	d7	d9	A	d10	d8	d6	d4	d2	ic	ic	ic	F	K	—	<b>B19A</b>
<b>9810KB 9813KB 9816KB</b>	ic	d1	d3	d5	d7	d9	d11	d13	ic	A	d14	d12	d10	d8	d6	d4	d2	ic	F	K	<b>B20-102</b>
<b>9811KB 9814KB 9817KB</b>	ic	d1	d1	d3	d5	d7	d9	d11	ic	A	d12	d10	d8	d6	d4	d2	ic	ic	F	K	<b>B20-102</b>
<b>9812KB 9815KB 9818KB</b>	ic	d1	d1	d1	d3	d5	d7	d9	ic	A	d10	d8	d6	d4	d2	ic	ic	ic	F	K	<b>B20-102</b>


# Magnetic Focusing Assemblies

## Types C121, C122

The EMI magnetic focusing assemblies are designed to provide a reduction of the effective cathode area and a corresponding decrease in the cathode dark current in end window photomultiplier tubes. This results in improvements in the signal/dark current ratio and so in the limit of detection of low intensity light. The two versions available, type C121 and C122, produce effective cathode diameters of approximately 6 mm and 12 mm respectively and are for use with 50 mm diameter photomultiplier tubes.

The magnetic field produced by the assembly deflects electrons from the edge of the cathode and the side walls away from the dynode system so that they are not multiplied. As a result, a considerable reduction in dark current is achieved.

The magnets can be used successfully where a small diameter incident beam is convenient or where the incident radiation can be focused on to a small area of the photocathode. They are designed for use with thermoelectric and dry ice photomultiplier tube housings manufactured by Products for Research Inc, as well as for custom built assemblies.



### CHARACTERISTICS

		Type C121	Type C122
<b>Outside diameter</b>	max. mm	51.8	51.8
<b>Depth of assembly</b>	max. mm	12.0	17.0
<b>Effective cathode dia. (see note 1)</b>	nom. mm	6.0	12.0
<b>Approx. dark current as a percentage of unfocused dark current (see note 5)</b>		5	20

### NOTES

- The nominal effective cathode diameter (see below) in millimetres is denoted by the figure after the stroke in the type number. The standard types are the C121/6 and the C122/12 providing effective cathode diameters of 6 mm and 12 mm respectively. The effective cathode diameter is the full width of half-maximum amplitude of the anode output as a small spot of light is swept across the cathode.
- The assembly must be positioned so that the surface which is flush with outer collar is against the photomultiplier window. It may then be held in place with adhesive tape. Alternatively, when the thermoelectric or dry ice housings, manufactured by Products for Research Inc, are to be used, the assembly should be fitted inside the corrugated pressure sleeve in the housing immediately in front of the photomultiplier window.
- A mu-metal shield should be used with the magnet for optimum performance and ideally both should be maintained at cathode potential. If this is not possible, the assembly should be left unconnected and should not be earthed if the anode is at earth potential.
- Care should be exercised when handling these magnets and they should not be stored in close contact with one another.
- The figures quoted for the 'focused' dark current as a percentage of the 'unfocused' dark current are for a 9558B operated at an overall sensitivity of 200A/lm. The percentage reduction obtained will vary from tube to tube and will depend on the particular type of photocathode employed.

# Information on other new tubes and other products

## Other new tubes

As a result of a continuous development programme and a policy to design tubes to suit the ever increasing number of applications, new variants of existing types may be added to the range at any time. At the time of printing this supplement, several new versions are planned for immediate introduction and it is hoped that the following brief details will serve as a guide until full information is printed. The parent types, where referred to, are described in the EMI photomultiplier catalogue reference P001 /fP70 available free on request.

- 9836** This is a 50mm (2 in) diameter tube having the same basic construction and specification as the 9656 series except that the end window is made from lime soda glass instead of borosilicate. It is intended for use in equipment where the lower boron content and thus lower neutron absorption of soda glass may be critical.
- 9769** This is a 50mm (2 in) diameter tube based on the 9656 but having 9 dynodes instead of 10. The main application is in scintillation counter assemblies where good linearity is required and high gain is not necessary. The tube will operate satisfactorily when plugged directly into sockets wired for the 9656 or other 10 stage tubes.
- 9840** This is a 50mm (2 in) diameter tube having dimensions and specifications similar to the 9750 but it incorporates an internal end window surface of triangular prisms as used in the 9778 described in this supplement. Like the 9778, the 9840 is designed for high resolution gamma cameras.

## Other products

### PHOTOMULTIPLIER TUBE POWER SUPPLY PM25A

A photomultiplier tube stabilised power supply will be available from EMI when this supplement is printed. It incorporates design features which we consider important for successful photomultiplier tube operation particularly under laboratory conditions. The following is an extract from the specifications :- output variable over the range 100-2500 volts at 5 mA — dual polarity, resettable to within 1 volt, output stability 1 part in  $10^5$  for 10% mains change, dimensions  $8\frac{1}{2}$  in. x 5 in. x 9 in., mains input 100/125 or 200/240 volts 48—66 Hz. For further information telephone Extension 2126 or write.

## Electron tube range

The EMI ELECTRON TUBE DIVISION manufactures a wide range of special electron tubes for equipment used in broadcasting, radar, nuclear and scientific applications.

PHOTOMULTIPLIER TUBES Ext. 2074  
Photomultiplier tubes which convert very low levels of illumination into usable electric currents are used extensively in astronomy, spectrophotometry, scintillation counting, spectrometry and broadcast television.

PHOTOMULTIPLIER TUBE HOUSINGS Ext. 2283  
A range of cooled and uncooled photomultiplier tube housings, including thermoelectric, dry ice and liquid nitrogen versions are available for optimum photomultiplier tube operation.

CAMERA TUBES Ext. 2078  
There is a wide range of vidicons, including all-electrostatic, available in various grades from general surveillance to broadcast studio.

IMAGE INTENSIFIERS Ext. 2075  
The image intensifier tube, capable of multiplying light up to a million times, is important for such applications as microscopy and astronomy.

CATHODE RAY TUBES Ext. 2073  
EMI activities in pioneering television have generated a range of specialised cathode ray tubes for radar and telecine work.

SPECIAL PRODUCTS Ext. 2551  
EMI manufactures the Printicon, a small all electrostatic monoscope; the Ebitron, a low light level intensifier-vidicon camera tube and spectroscopic lamps. Two types of spectroscopic lamp are available, hollow cathode and electrodeless discharge tubes together with a microwave power generator. A range of printed circuit scanning coils and complete scanning assemblies for 13 mm, 26 mm and 30 mm vidicon camera tubes is also produced.

SOLID STATE PHOTODIODES Ext. 2126  
These include a range of linear and avalanche silicon photodiodes including fast and rugged types having wide spectral response.

PRECISION MICROMESH Ext. 2073  
The very fine metallic mesh currently employed in EMI vacuum tubes is also used in various other branches of industry and science, such as microscopy, mass spectrometry, biology, filtering and optics.

The EMI Electron Tube Division has great experience and comprehensive facilities in research, development and manufacture of light sensing and light emitting devices and allied equipment.



**EMI Electronics Ltd Electron Tube Division**  
243 Blyth Road, HAYES, Middlesex, UB3 1HJ, England.  
Telephone: 01-573 3888 Extension 2542  
Cables: Emidata, London Telex: London 22417

